

Assessment and Comparison of LAMP and PCR Techniques in Diagnosis of *Cryptococcus neoformans*

Mohammad Hassan Shahhosseiny^{1*}, Mohammad Ghahri² and Elham Moslemi³

¹Department of Microbiology, Islamic Azad University, Qods branch, Tehran, Iran.

²Department of Biology, Imam Hossein University, Tehran, Iran.

³Iranian Gene Fanavar institute (IGF), Tehran, Iran.

(Received: 11 July 2013; accepted: 24 September 2013)

The purpose of this study is comparative assessment of traditional methods and two molecular techniques in diagnosis of *C. neoformans* in cerebrospinal fluid (CSF) as well as culturing samples. 300 CSF samples from patients with suspected meningitis have been collected from multiple hospitals. Samples were evaluated by staining, culture and PCR and LAMP. The target gene of primers used in this study for PCR, selected on the basis of 16S rRNA, also, the LAMP test primers were designed with the help of "Primer explorer V4" software based on the *ura5* gene. Internal Control of PCR constructed by the competitive strategies and cloned in pTZ57R. The PCR results were positive in 4 cases, while the results of culturing and staining tests were negative for all the CSF samples. The sensitivity of LAMP test was determined up to 5 cells of *C. neoformans*. The LAMP test was positive in 5 cases of CSF samples (four positive cases of PCR and another sample with a negative PCR). The LAMP technique which provides additional advantages such as higher sensitivity and specificity and no need to expensive hardware, also eliminates the identification stage by electrophoresis methods, and available especially in developing countries.

Key words: PCR, LAMP, *Cryptococcus neoformans*.

C. neoformans is a yeast member of Sporidiobolaceae family with a polysaccharide capsule which often replicates by budding¹. Based on the capsular polysaccharide antigen, the fungi are categorized into five serotypes, including A, D, AD, B and C. This species composed of two varieties according to the molecular, biochemical, ecological and epidemiological differences namely: *C. neoformans* var. *neoformans* which includes A, D and AD serotypes and *C. neoformans* var. *gattii* which includes B and C serotypes^{2,3}. Recently,

scientists called the A serotype and the D serotype as *C. neoformans* var. *grubii* and *C. neoformans* var. *neoformans*, respectively⁴⁻⁵. Studies suggest a close similarity between the *C. neoformans* var. *gattii* and the *C. neoformans* var. *neoformans* which composed of two non-pathogenic species namely *C. amyloletus* and *Tsuchiyaea wingfieldii* as well as *Filobasidiella depauperata*⁶.

In the absence of strong immune system capability, Cryptococcosis acts as a human opportunistic infectant, furthermore in current decades, regarding to increasing of HIV positive people it has become very important in public health level⁷. Using more sensitive and specific laboratory techniques for rapid detection and early antifungal treatment in diagnosis of neurological cryptococcosis, is essential (8).

* To whom all correspondence should be addressed.
Tel: +98 21-44844946; Fax: +98 21-44861889;
E-mail: shahhosseiny@yahoo.com

C. neoformans infections spread out through inhalation of the dust contaminated with living fungal conidia in the human, and usually people with cellular immune system deficiency will be infected⁹. The common appearances of disease are as lung and cervical forms, but due to spreading out of fungi through the primary pulmonary centre, other various forms of disease were reported as cutaneous, mucosal, bone and visceral. The human pathogen species is *C. neoformans*, which prefers to attack the central nervous system, because of unknown reasons. The infection is spread out to these centers, by apparent or unapparent pulmonary infections¹⁰.

Generally, the culture-based methods are time-consuming techniques with low sensitivity which require both adequate facilities and experienced specialists to analyze the results; in addition, there is not enough sensitivity in staining methods. Compared to aforementioned two methods, serological methods which are based on antigen - antibody reaction, as well as using monoclonal antibodies, and the flow cytometry methods have higher sensitivity and accuracy¹¹. However, false positive and/or negative results could be achieved due to cross reactions. Thus, we need a method with high sensitivity; specificity and speed performance to identify this agent needed which has no any limitations of above methods.

Previous studies suggested; because of their higher sensitivity, specificity and precision, the diagnostic methods based on PCR and the loop mediated isothermal amplification (LAMP) test sound to be appropriate methods to detect pathogens such as *C. neoformans* today^{8,12}. The level of PCR sensitivity and accuracy in some specific parts of *C. neoformans* genome recognition depends on several variables such as DNA extraction method, the target gene, the designed primers and the method of product identification which are different in various studies¹³.

The prepared primers are generally designed for following common portions: 16S rRNA gene¹⁴, capsular gene¹⁵ or other targets such as phospholipase enzyme gene¹⁶. However, to design the species-specified primers, more emphasis was on the existing constant parts in the rRNA gene sequences of *C. neoformans*. Despite of positive features of PCR method, some

disadvantages such as expensive equipments and other problems make its widespread usage difficult, especially in routine laboratories.

Because of various reasons PCR method could be replaced, especially in developing countries, with Loop-mediated isothermal replication technique (LAMP), as a suitable method in diagnostic purposes¹⁷. This technique (LAMP) is a simple and time saving method which totally identifies 6 or 8 regions in the DNA model. Typically, compared to the other existing DNA isothermal replication methods, this test has very high sensitivity and specificity^{18,19,20}, so that large amount of DNA (10 -30 micrograms in 25 µl) can be made in a short period of time (30-60 minutes) maintaining high specificity²¹. Nowadays, the LAMP method has been commercially used to identify various pathogens. It can be considered as a significant technique in developing countries¹⁷.

Several diagnostic methods compared in many studies which have been conducted in the diagnosis of *C. neoformans*. For example, three cultural, serologic and PCR methods compared by Rappelli et al . According to the results, the sensitivity of cultural, serologic and PCR methods in diagnosis of *C. neoformans* are 67%, 89% and 98%, respectively²². Rezai et al 2012, to detect the fungi in the serum of HIV positive people, used LAMP method and reported four positive cases.

This study aims to assess two molecular techniques; PCR and LAMP as two very practical methods in diagnose of *C. neoformans* in the infected samples such as cerebrospinal fluid as well as culture samples and comparing them.

MATERIALS AND METHODS

C. neoformans culturing: The provided strain obtained from the Microbiological Collection Center of Scientific and Industrial Research Organization (ATCC number: 13690) and two provided strains from the Resalat laboratory, Tehran (isolated from patient) were cultured on the Niger seed agar medium (Acumedia Company) and the plates were incubated for a week at 37 ° C (*C. neoformans* colonies usually grow within three days). DNA was extracted from three strains and the PCR and LAMP tests were optimized on them.

Sample collection

300 samples of cerebrospinal fluid (at least 500 μ l) from patient suspected to meningitis have been collected during 2012-2013, which then transferred to the laboratory under cold conditions. The samples were centrifuged at 8000 rpm for 2 minutes. Then, the supernatant fluid was slowly evacuated; and remaining sediment was solved in 120 μ l PBS buffer and 40 μ l of the prepared solution was placed on the lamel and mixed with a drop of Hindi ink. The sample was observed under microscope after placing a lamella on it, in case of existence of *C. neoformans* a transparent and clear zone can be seen around the yeast cells which appeared as a sky scene full of stars. Another 40 μ l of the sediment were cultured on the Niger seed agar medium (Acumedia Company) as well, and the plates were incubated at 37 °C for a week (*C. neoformans* colonies usually grow within three days).

DNA extraction from cerebrospinal fluid samples

To extract the DNA, DNP Cinaclon kit (Cat: DN811540) was used according to the kit instructions. The genomic DNA of bacteria, yeast and humans used in this study for the characteristics tests was also extracted using same method.

Primer design for PCR and LAMP

The used primers in this study for PCR technique are 16S rRNA gene-specific primers of *C. neoformans* (Table 1)²³. Also, the LAMP test primers were designed with the help of "Primer explorer V4" software based on the existing sequence in the *ura5* gene (ACCESSION EU399581)(primer explorer V4; <http://primerexplorer.jp/E/>) (Table 1).

PCR reaction

5 μ l of the extracted DNA using the ITS1 and CN4 primers specific to 16S rRNA gene with the final concentration of 0.4 μ l (1 μ l of 10 mM concentration) were duplicated in the final volume of 25 μ l by mixing with the followings: 2.5 μ l of PCR (10X) (Cinaclon) buffer; 0.75 μ l of magnesium chloride ($MgCl_2$) (of 50 mM concentration (Cinaclon)); 0.5 μ l of dNTP mixture (10 mM) (Cinaclon) and 0.3 μ l of Taq DNA Polymerase enzyme (Cinaclon). The thermal cycles consisted of a six-minute cycle at 94 °C, 40 sequential cycles (including 30 seconds at 94 °C, 30 seconds at 62 °C and 1 minute at 72 °C,) and a final polymerization

cycle of 72 °C for 10 minutes. After ending the thermal cycles, the electrophoresis technique was used to observe the amplified DNA. The PCR product alongside the size marker and positive and negative controls was studied on 1.5% agarose gel by using ethidium bromide or SYBR Green I and UV light on the Transilluminator device.

Construction of internal control (IC) and Ideal density

To construct the competitive internal Control of PCR test, we re-designed and synthesized the forward and reverse primers of Cryptococcus PCR at 5' portion of the kinetoplast gene primers of *Leishmania* with the product size of 620 pairs (Table 1). Doing the PCR reaction to replicate the internal control fragment, 5 μ l of *leishmania*'s extracted DNA the front and rear IC primers with the final concentration of 0.4 mM were amplified in the final volume of 25 μ l were used by mixing with the followings: 2.5 μ l of 10X PCR buffer; 0.75 μ l of magnesium chloride ($MgCl_2$ 50 mM); 0.5 μ l of mixed dNTP (10 mM) and 0.3 μ l of Taq DNA Polymerase enzyme (5 unit/). The thermal cycles consisted of 35 sequential cycles (including 30 seconds at 94 °C, 30 seconds at 60 °C and 1 minute at 72 °C,) and a final polymerization cycle of 72 °C for 25 minutes. After ending the thermal cycles, the electrophoresis technique was used to observe the amplified DNA. By using Blue/White Screening the resulting construction was selected. The suitable clones finally extracted via plasmid extraction and PCR test.

Different concentrations of standard *C. neoformans* IC and DNA were examined In order to obtain the optimal density of the plasmid containing IC which not competing the target DNA, The DNA concentration of *C. neoformans* as well as the IC concentration used in the test was determined by a spectrophotometer at A260. For PCR test, 5 μ l of the DNA obtained from the samples were mixed with 1 μ l of IC, and was used in the PCR mixture.

LAMP reaction

The reaction was performed at the volume of 25 μ l, including the following reagents: a mixture of 5.2 μ l twice-distilled water, 2.5 μ l of enzyme buffer with the 10X concentration, 1 μ l of a mixture of external primers (F3/B3) with a concentration of 0.2 μ M and internal primers (BIP/FIP) with a concentration of 1.6 μ M, 1 μ l of loop

primers (LB/LP) with a concentration of 0.8 µM, 1.8 µl of MgSO₄ with a concentration of 7 mM, 1 µl (8U) of Bst DNA polymerase enzyme (New England biolab), 4 µl of Betain (Sigma) with a concentration of 0.8 M, 3.5 µl of dNTP with a concentration of 1.4mM and 5µl of the extracted DNA. The reaction was performed at the temperature of 66° C within 60 minutes. Positive and negative controls were used in each round of the reaction. To evaluate the reaction product, 1 µl of 10 times diluted SYBR green (Invitrogen lot: 49753A) was added to each reaction tube and was reviewed on the Transilluminator device at the 302 nm wavelength. The tubes containing positive and negative reactions were observed having fluorescent green and orange colors, respectively.

Sensitivity and specificity of LAMP and PCR reaction

To determine the sensitivity of optimized LAMP and PCR tests, different dilutions of the *C. neoformans* in deionized water were prepared in the light of direct counting microscopic method using extracted DNA from different dilutions. Further, to determine the specificity of the optimized PCR and LAMP tests, their functions assessed on the DNAs of human, mice, *Candida albicans*, *Aspergillus niger*, *Escherichia coli*, *Fusarium solani*, *Saccharomyces cerevisiae*, *Pichia pastoris* and a number of other agents, along with positive and negative control samples.

PCR product and internal control Cloning

After PCR product and IC purification, both amplicons were cloned at pTZ57R vector using the T/A cloning Thermo scientific kit(K1214). The resulting plasmids were extracted by the Plasmid Mini Extraction Kit of Bioneer

Company(K3112). Then, the plasmids containing the PCR product and the IC were confirmed, using PCR method.

Sequence finding

Using “Dideoxy Chain Termination” method by a Macrogen inc, DNA sequencing was conducted in two directions

RESULTS AND DISCUSSION

The PCR test was optimized on DNA extracted from the *C. neoformans* strain (cultured on the medium) (Figure1A). The diagnosis limitation (sensitivity) was obtained at least as 10 copies in the examined sample by preparation series of dilutions of the yeast culture and performing DNA extraction and PCR test (Figure 1B). Also, no unwanted product was seen in the specificity test.

In Figure 2, the results of internal control PCR test and its examination by different DNA amounts of *C. neoformans* are seen which have variety of 100 to 3 million fungi. The IC ideal number for PCR test was obtained as 1000 per reaction and then used for disease samples examination.

To diagnosis of Cryptococcus infections in the cerebrospinal fluid samples, specific ITS1 and CN4 primers were used. While the results of culturing and staining tests were negative for all the CSF samples PCR test results were positive in 4 cases of 300 CSF samples. To review and approve the replicated PCR product in the disease sample, the 415-bp fragment was sequenced by “Sequencing” method. As well, nearly 100 percent compliance appeared by comparing the resulting sequence with existing sequences through BLAST application.

Table 1. Used primers in PCR and LAMP tests and making internal control

Primer	Sequence (5'.....3')
ITS1	5' TCC GTA GGT GAA CCT GCG G 3'
CN4	5' ATC ACC TTC CCA CTA ACA CAT T 3'
ICF	5' TCC GTA GGT GAA CCT GCG G -TCG CAG AAC GCC CCT ACC 3'
ICR	5' ATC ACC TTC CCA CTA ACA CAT T -AGG GGT TGG TGT AAA ATA GGC 3'
F3	5' TCC TTG GCT GCT GTC TCC 3'
B3	5' GCC TTG CCA GAG GTA AGA AC 3'
FIP	TGG GAC AGA CTC ACG TCC TTC TAA ACC GGC AAA GAT ATC GGC 3'5'
BIP	5' CCA GTG CGA CAG CGA TGA GCT TCG TCC CTT GAG AGG CG 3'
LF	5' TCTCCTTCCTGTTGTAGCAGTA 3'
LB	5' AGCCAGTAGCACGGTGAGGG 3'

The LAMP reaction was optimized on extracted DNA from the standard strain and two isolated strains from patients at a temperature of 66 °C within 1 hour (Figure 3). The sensitivity of LAMP test with optimized parameters was determined up to 5 cells of *C. neoformans* (Figure 4). Also, very high specificity appeared in the LAMP test, since it acts specifically using six specific primers so that it only reacted with *C.*

neoformans DNA and no reactions occurred with DNA of other used agents (Figure 5).

In 5 cases of CSF samples, the LAMP test was positive (four positive cases of PCR and another sample with a negative PCR).

Early and accurate diagnosis of infectious agents are Factors which reduce mortality, morbidity and the costs of hospital infectious diseases as well. Although the traditional




Fig. 1(a). The optimized PCR test on the standard *Cryptococcus neoformans* strain

M line: Fermentas size marker (100 bp DNA Ladder); Line 1: Positive PCR test;

Line 2: Negative control;

Fig. 1(b). PCR sensitivity test using counted colonies of *Cryptococcus neoformans* (indicated CFU)

M line: Fermentas size marker (100 bp DNA Ladder); Line 1: Positive control samples;

Lines 2-6: Contain certain samples, including, respectively: 2: (CFU 100000), 3: (CFU 10000), 4: (CFU 1000), 5: (CFU 100), 6: (CFU 10), 7: (CFU 5), 8: (CFU 1), 9: Negative control sample




Fig. 2. The PCR electrophoresis results with certain number of *Cryptococcus neoformans* DNA:

1: Fermentas 100 bp size marker; 2: Positive control (containing IC and *Cryptococcus* DNA); 3: DNA of a number of 3000000 agents; 4: A number of 2000000; 5: A number of 1000000; 6: A number of 100000; 7: A number of 10000; 8: A number of 1000; 9: A number of 100; 10: Negative control (merely containing IC).




Fig. 3. The optimized LAMP test.

Tubes 1-3: Positive reaction with three strains of *Cryptococcus neoformans*;

Tube 4: Negative reaction

culturing-based methods are still considered as a gold standard in some cases, but in case of need to rapid identification of microorganisms, including *C. neoformans* they generally doesn't satisfy characteristics of an optimal method. False negative results in culture and staining methods, which need living cells, is due to their low

sensitivity and sometimes the samples may be received late to the laboratory, so cells may have been destroyed. Time detection of these yeasts in patients with neuro-cryptococcosis plays an important role in reduced mortality rate from this disease. In Iran and other developing countries the low rate isolation of *C. neoformans*, compared




Fig. 4. Determining the sensitivity of LAMP test

4a. Tubes 1: positive control; 2: 1/000/000 cells; 3: 100/000 cells; 4: 10/000 cells; 5:1000 cells; 6: 100 cells; 7: 50 cells; 8: 5 cells; 9: Negative control.

4b. Lines 1: 100 bp DNA Ladder Plus (Fermentas) size marker; 2: Positive control; 3: 1/000/000 cells; 4: 100/000 cells; 5: 10/000 cells; 6:1000 cells; 7: 100 cells; 8: 50 cells; 9: 5 cells; 10: Negative control.




Fig. 5. Determining the specificity of LAMP test

Tube No. 1: positive control; Tube No. 2: LAMP reaction with human DNA; Tube No. 3: With *Candida albicans* DNA; Tube No. 4: With *E. coli* DNA; Tube No. 5: With *Saccharomyces cerevisiae* DNA; Tube No. 6: With *Aspergillus niger* DNA; Tube No. 7: With *Fusarium solani* DNA; Tube No. 8: With *Pichia pastoris* DNA; Tube No. 9: Negative control

to other countries does not indicate low prevalence of infections caused by this yeast, but the difference is related to identification and diagnosis methods. However, studies have shown higher proportion of *C. neoformans* will be detected using much more accurate and sensitive detection methods.

To identify *C. neoformans*, several methods have been already suggested, each one has some advantages and disadvantages. Studies using cultivation techniques along with molecular methods such as the PCR and LAMP methods have shown the culturing method is time-consuming and appears high false negative results. Also, cultures' result will be reported as negative in cases which treatment is done before sampling. So, culturing method results insufficient data in *C. neoformans* detection and application of molecular methods considered as the first priority in diagnosis the agent⁸. Further, more sensitivity and specificity of molecular methods clearly proved by several studies compared serological with molecular methods, such as PCR to detect *C. neoformans*,²⁵.

In 1984, Cohen, reported that occurrence percentage of false negative results is more than false positive results by comparing diagnostic methods of "Indian Ink" staining, "Acridine Orange" fluorescent staining and serologic Latex agglutination method on 162 *C. neoformans* suspected samples²⁶. Velegraki et al. 2001, concluded PCR test on *ura5* gene followed by performing the RFLP technique using an "In House" or commercial DNA extraction method can be a good strategy to fast detect of different varieties of cryptococcus by using PCR technique and applying specific primers of *C. neoformans ura5* gene,²⁷.

In 2009, Klein et al. reviewed Glycine consuming test results on the CGB Agar medium (l-canavanine glycine bromothymol blue agar) and urea hydrolysis test and melanin production on 102 isolated samples of *C. neoformans* var. *gattii* simultaneously and compared them to results obtained by the "D2 LSU DNA Sequencing" method (by using the "rDNA" of ribosomal large subunit). They declared application of CGB Agar medium with urea hydrolysis test and melanin production on Niger seed Agar can provide good results to distinguish two mentioned species just

like as the "D2 LSU DNA Sequencing" method²⁸.

To extract this agent DNA, various methods have been used in different studies including commercial kits such as QIAGEN²⁹, Breeden Lab²³, using lysis buffer³⁰, and freezing / boiling method³¹.

There is a little Study which compares different methods of DNA extraction of *C. neoformans* in molecular diagnostics. In present study, the DNP kit was used for clinical samples obtaining good results.

Various parts of the *C. neoformans* genome such as 16S rRNA gene¹⁴, polysaccharide capsule gene¹⁵ and phospholipase enzyme gene¹⁶ incorporated in molecular methods such as PCR under different protocols and techniques to detect the agent. The 16S rRNA gene of genome used by PCR and the *ura5* gene were used in the LAMP technique as well because of specificity and constant sequence of that region to detect the agent. In molecular typing of *C. neoformans* and identifying their varieties the followings were used (24).

Finally, sufficiency of ITS1 and CN4 primers to detect *C. neoformans* in cerebrospinal fluid samples as well as, the designed LAMP primers were proved. The optimized conditions of this study led to produce a specific product with high degree of specificity and sensitivity and no unwanted by-products. Although the number of positive samples in this study was not significant, but the existence of cryptococcus or at least its genome in the tested CSF samples appeared in PCR and LAMP tests, which the results were negative in direct microscopic and culturing experiments. The study results also confirmed this issue is true. Although only four positive samples in PCR test and five positive samples in the LAMP test were obtained, but we did not also succeed to identify *C. neoformans* in any of the CSF samples using the available culturing and staining methods. In another similar researches which staining, culturing and PCR methods conducted to detect the yeast, the number of positive samples in PCR were further than culturing and staining methods, which indicates more sensitivity and accuracy of PCR comparing the culturing and staining methods (32).

This research suggests LAMP technique is sufficient. First time in 2000, the LAMP technique

was introduced by a Japanese scientist, Notomi, and his colleagues. Since this technique is able to detect the target sequence in the 6 separated regions by 4 specific primers of FIP, BIP, F3, B3, it acts so specifically in target sequencing. Also, there are more advantages of this method compared to other molecular techniques such as PCR, NASBA, 3SR and SDA (19). Nagamin in 2002, added two primers called loop primers (FLP and BLP) to the reaction and thus the synthesis of DNA in this reaction became faster (18). In 2004, the Japanese scientist, Ushikubo, stated that adding two loop primers of FLP and BLP to the LAMP technique can reduce the required time of the test from one hour to less than half an hour, and also, a great amount of white precipitate of pyrophosphate magnesium is produced as the accessory product in the experiment, which even makes the diagnosis of a positive test possible visually (34). Isothermal replication using a loop is a particular gene duplication method, which DNA is duplicated just by an enzyme with DNA strand displacement feature (such as Bst polymerase) under isothermal conditions (19,33). The frontal inner primer (FIP) includes F2 fragments and the F1 complementary sequence (F1c), and the back inner primer (BIP) includes B2 fragments and the B1 complementary sequence (B1c) in this study, which sometimes allocates larger space sequences of thymidine between these dual fragments that are not necessary either (21,35); F3 and B3 are paired with the initial pattern strand, and with the help of enzyme property to pull the strand, a two-strand is completed and extracted, and simultaneously the strand produced by F1c and B1c will remain as a single-strand for the successive cycles of F1c, B1c, LB, LF performance (19,42). Simply, we could see the LAMP test result by the naked eye. So, this observation is obtainable by using two methods, whether through creating turbidity by addition of cationic polymers or through the addition of SYBR Green fluorescent dye to the LAMP products (20,21). It was proved in this study to detect the opportunistic yeast in clinical samples such as CSFs suspected to neurocryptococcosis, the LAMP technique is more sensitive than the PCR technique. Meanwhile, low cost of this technique regarding required hardware and materials, introduced it as a good alternative for PCR. In several studies other agents such as B and C

hepatitis (36, 37), *Mycobacterium tuberculosis* (38) and many other agents were detected, and capabilities of the LAMP technique comparing other molecular methods, was proved, especially in diagnostic purposes.

The study results clearly indicate PCR assay based on 16S rRNA gene sequences and the LAMP test based on *ura5* gene are both valuable and reliable techniques for *C. neoformans* detection in cerebrospinal fluid. However, the LAMP technique is a more available method with additional advantages such as higher sensitivity and specificity has no need to expensive hardware and eliminating the identification stage by electrophoresis methods, especially in developing countries.

ACKNOWLEDGEMENTS

Hereby we would like to appreciate the Islamic Azad University- shahr-e-Qods branch for providing the expenses of this study. Special thanks from Iranian Gen Fanavar institute (IGF) particularly Dr. Parastoo Chamanrokh, Sam Masahi, Amin Mahmudi, Atefeh Shahhosseiny for their scientific generosity and providing laboratory equipments.

REFERENCES

1. Idnurm, A., Bahn, Y.S., Nielsen, K., Lin, X., Fraser, J.A., Heitman, J. Deciphering the Model Pathogenic Fungus *Cryptococcus Neoformans*. *Nature Reviews Microbiology* .2005;**3**, 753-764.
2. Kwon-Chung, K.J., Varma, A. Do major species concepts support one, two or more species within *Cryptococcus neoformans*? *FEMS Yeast Res* 2006; **6**: 574–587.
3. Okabayashi, K., Kano, R., Nakamura, Y., Watanabe, S., Hasegawa, A. Capsule-associated genes of serotypes of *Cryptococcus neoformans*, especially serotype AD. *Med Mycol* 2006; **44**: 127–132.
4. Nielsen, K.DOA., Heitman, J. *Cryptococcus neoformans* mates on pigeon guano: implications for the realized ecological niche and globalization. *Eukaryot Cell*. 2007; **6**:59-949.
5. Bovers, M., Hagen, F., Boekhout, T. Diversity of the *Cryptococcus neoformans*- *Cryptococcus gattii* species complex. *Revista Iberoamericana de Micologia.*, 2008; **25**: S4–12.
6. Findley, K.R-CM., Metin, B., Kroiss, J., Fonseca, A., Vilgalys, R., Heitman, J. Phylogeny

- and phenotypic characterization of pathogenic *Cryptococcus* species and closely related saprobic taxa in the Tremellales. *Eukaryot Cell*. 2009; **8**(3): 61-353.
7. Pan, W., Khayhan, K., Hagen, F., Wahyuningsih, R., Chakrabarti, A., Chowdhary, A., Ikeda, R., Taj-Aldeen, S.J., Khan, Z., Imran, D., Sjam, R., Sriburee, P., Liao, W., Chaicumpar, K., Ingviya, N., Mouton, J.W., Curfs-Breuker, I., Boekhout, T., Meis, J.F., Klaassen, C.H.W. Resistance of Asian *Cryptococcus neoformans* Serotype A Is Confined to Few Microsatellite Genotypes. *PLoS One.*; **7**(3): e32868.
 8. Rezaei-Amirabadi, A., Shahhosseiny, M.H., Vand-Yousefi, J., Ghahri, M. Loop mediated isothermal amplification (LAMP) in Diagnosis of neurocryptococcosis. *African Journal of Biotechnology*. 2012; **11**(17): 3986-3992.
 9. Chen, J., Varma, A., Diaz, M.R., Litvintseva, A.P., Wollenberg, K.K., Kwon-Chung, K.J. *Cryptococcus neoformans* Strains and Infection in Apparently Immunocompetent Patients, China. *Emerging Infectious Diseases*. 2008; **14**(5):755-762.
 10. Kwon-Chung, K.J., Varma, A. Do major species concepts support one, two or more species within *Cryptococcus neoformans*? *FEMS Yeast Res* 2006; **6**: 574–587.
 11. Cleare, W., Casadevall, A. The Different Binding Patterns of two Immunoglobulin M Monoclonal Antibodies to *Cryptococcus neoformans* Serotype A and D Strains Correlate with Serotype Classification and Differences in Functional Assays. *Clin Diagn Lab Immunol*. 1998; **5**(2): 125–129.
 12. Raimondi, A., Ticozzi, R., Salia, G., Grazia-Blotti, M. Genotype-based differentiation of the *Cryptococcus neoformans* serotypes by combined PCR-RFLP analysis of the capsule-associated genes CAP10 and CAP 59. *Medical Mycology*. 2007; **45**: 491-501.
 13. McFadden, D.C., Casadevall, A. Capsule and melanin synthesis in *Cryptococcus neoformans*. *Med Mycol*. 2001; **39** Suppl 1:19-30.
 14. Pashoal, R.C., Hirata, M. Neurocryptococcosis: Diagnosis By PCR Method. *Revar. Inst. Med. trop. S. Paulo*. 2004; **46**(4): 203-207.
 15. Chang, Y.C., Kwon-Chung, KJ. Isolation, characterization, and localization of a capsule-associated gene, CAP10, of *Cryptococcus neoformans*. *J Bacteriol*. 1999; **181**: 5636-5643.
 16. Anastasia, P., Litvintseva, L., Kestenaum, R. Comparative Analysis of Environmental and Clinical Population of *Cryptococcus neoformans*. *Journal of Clin Microbiol*. 2005; **43**(2): 556-564.
 17. Mori, Y., Notomi, T. Loop-mediated isothermal amplification (LAMP): a rapid, accurate, and cost-effective diagnostic method for infectious diseases. *J Infect Chemother*. 2009; **15**:62–69.
 18. Nagamine, K. Accelerated reaction by loop-mediated isothermal amplification using loop primers. *Molecular and Cellular Probes*. 2002; **16**(3): 223-229.
 19. Notomi, T., Okayama, H., Masubuchi, H., Yonekawa, T., Watanabe, K., Amino, N. Loop-mediated isothermal amplification of DNA. *Nucleic Acids Res*. 2000; **15**: 28(12):E.63
 20. Parida, M.S.S., Kumar-Dash, P., Rao, P.V.L., Morita, K. Loop mediated isothermal amplification (LAMP): a new generation of innovative gene amplification technique; perspectives in clinical diagnosis of infectious diseases. *Rev Med Virol*. 2008; **18**: 407-21.
 21. Mori, Y., Hirano, T., Notomi, T. Sequence specific visual detection of LAMP reactions by addition of cationic polymers. *BMC Biotechnol*. 2006; **6**:3.
 22. Rappelli, P., Are, R., Casu, G., Fiori, PL., Cappuccinelli, P., Aceti, A. Development of a nested PCR for detection of *Cryptococcus neoformans* in cerebrospinal fluid. *J.Clin. Microbiol*. 1998; **36**: 3438-3440.
 23. Kwon-Chung, KJ., Edman, J.C., Wickes, B.L. Genetic association of mating types and virulence in *Cryptococcus neoformans*. *Infect. Immun*. 1992; **60**: 602-605.
 24. Meyer, W.K.S., Castañeda, A., Jackson, S., Huynh, M., Latouche, G.N., editor. Global molecular epidemiology offers hints towards ongoing speciation within *Cryptococcus neoformans*. In: Abstracts of the 5th International Conference on *Cryptococcus* and Cryptococcosis; 2002 March 7-3,2002; South Australian Postgraduate Medical Education Association, Adelaide, Australia.
 25. Dolan, C.S., Immaculata, X., Ashutosh, B.I., Dipankar, M., Bhowmik, V.A.R. Detection of *Cryptococcus* by conventional, serological and molecular methods. *J Med Microbiol*. 2009; **58**: 1098-1105.
 26. Cohen, J. Comparison of the sensitivity of three methods for the rapid identification of *Cryptococcus neoformans*. *J Clin Pathol*.1984; **37**(3):4-332.
 27. Velegraki, A., Kansouzidou, A., Smilakou, S., Mitroussia-Ziouva, A., Legakis, N.J. Prospective use of RFLP analysis on amplified *Cryptococcus neoformans* URA⁵ gene sequences for rapid identification of varieties and serotypes in clinical samples. *Medical Mycology*. .2001; **39**(5): 409-417.

28. Klein, K.R., Deml, S.M., Rysavy, J.M., Wohlfiel, S.L., Wengenack, N.L. Identification of *Cryptococcus gattii* by use of L-canavanine glycine bromothymol blue medium and DNA sequencing. *J Clin Microbiol.* 2009; **47**(11): 3669-72.
29. Mitchell, T.G., Perfect, J.R. Cryptococcosis in the era of AIDS—100 years after the discovery of *Cryptococcus neoformans*. *Clin. Microbiol. Revar.* 1995; **8**: 515-548.
30. Stephen, C., Lester, S., Black, W., Fyfe, M., Raverty, S. Multispecies outbreak of cryptococcosis on southern Vancouver Island, British Columbia. *Can. Vet. J.* 2002; **43**: 792-794.
31. Brandt, M.E., Hutwagner-Klug, L.C., Baughman, S., Rimland, D.E., Graviss, A., Hamill, R.J., Thomas, C., Pappas, P.G., Reingold, A.L. Pinner and the Cryptococcal Disease Active Surveillance Group. Molecular subtype distribution of *Cryptococcus neoformans* in four areas of the United States. *J. Clin. Microbiol.* 1996; **34**: 912-917.
32. Tanaka, K., Miyazaki, T., Maesaki, S. Detection of *Cryptococcus neoformans* gene in patients with pulmonary cryptococcosis. *J Clin Microbiol.* 1996; **34**: 2826–2828.
33. Nagamine, K., Hase, T., Notomi, T. Accelerated reaction by loop-mediated isothermal amplification using loop primers. *Molecular and Cellular Probes* .2002; **16**: 223-229.
34. Ushikubo, H. The principle of LAMP method. *Uirusu.* 2004; **12-54**:107.
35. Hong, T.C., Cuong, D.V., Parida, M., Minekawa, H., Notomi, T., Hasebe, F., Morita, K. Development and evaluation of a novel loop-mediated isothermal amplification method for rapid detection of severe acute respiratory syndrome coronavirus. *J Clin Microbiol.* 2004; **42**: 1956-61.
36. Moslemi, E., Shahhosseiny, M.H., Javadi, G.h., Praivar, K., Nejad-Sattari, T., Keyvani-Amini, H. Loop mediated isothermal amplification (LAMP) for rapid detection of HBV in Iran. *African Journal of Microbiology Research.* 2009; **3**(8): 439-445.
37. Nasr-Esfahani, S.h., Shahhosseiny, M.H., Yaghmai, P., Praivar, K., Moslemi, E. Keyvani-Amini, H. Rapid and simple detection of Hepatitis C virus by reverse transcriptase -loop-mediated isothermal amplification method. *African Journal of Microbiology Research* .2010; **4**(23): 2580-2586.
38. Kohan, L., Shahhosseiny, M.H., Razavi, M.R., Parivar, K., Moslemi, E., Werngren, J. Evaluation of loop mediated isothermal amplification for diagnosis of *Mycobacterium tuberculosis* complex in clinical samples. *African Journal of Biotechnology.* 2011; **10**(26): 5096-5101.